

Anti-WFDC2 antibody (50-100) {FITC} (STJ500963)

STJ500963

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-WFDC2 (50-100) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.68 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:750 DB: 1:20, 000 ELISA: 1:20, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10406
Gene Symbol	WFDC2
Uniprot ID	WFDC2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 50-100 on human WAP four-disulfide core domain protein 2.
Immunogen Region	50-100
Specificity	
Immunogen Sequence	