

Anti-PTK2 antibody (740-790) (STJ500981)

STJ500981

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Chicken polyclonal antibody anti-PTK2 (740-790) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Chicken
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.6 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10,000 IP: 1:200 IHC: 1:100 ICC: 1:100 IF: 1:100
Formulation	
Isotype	IgY
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5747
Gene Symbol	PTK2
Uniprot ID	FAK1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 740-790 on human focal adhesion kinase 1 protein.
Immunogen Region	740-790
Specificity	
Immunogen Sequence	



Immunohistochemistry of rat tissue with Anti-PTK2 primary antibody (STJ500981) at 1:100 dilution in Buffer. Section was treated with DAB and Haematoxylin stain. Magnification is at 5X.