

Anti-FBXW4 antibody (240-290) {Biotin} (STJ500991)

STJ500991

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-FBXW4 (240-290) is suitable for use in Dot Blot, ELISA, Immunoprecipitation and Western Blot research applications.
Applications	DB, ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.68-0.7 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:4, 000 ELISA: 1:4, 000 IP: 1:150
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	6468
Gene Symbol	FBXW4
Uniprot ID	FBXW4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 240-290 on human FBX4 protein.
Immunogen	240-290
Region	
Specificity	
Immunogen	
Sequence	