

Anti-FeSODB1 antibody (STJ501002)

STJ501002

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Chicken polyclonal antibody anti-FeSODB1 is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Chicken
Reactivity	L.infantum

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgY
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Purified peptide from FeSODB1 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	



Western blot with Anti-FeSODB1 Antibody (STJ501002) with recombinant protein. Antibody diluted at 1:500 in buffer. 500 ng and 1 µg protein/lane.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081