

Anti-FFAR2 antibody (250-300) (STJ501008)

STJ501008

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-FFAR2 (250-300) is suitable for use in Dot Blot, ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	DB, ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.58-0.59 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150 IHC: 1:50-1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	2867
Gene Symbol	FFAR2
Uniprot ID	FFAR2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 250-300 on human free fatty acid receptor 2 protein.
Immunogen Region	250-300
Specificity	
Immunogen Sequence	



Western blot with Anti-FFAR2 Antibody (STJ501008) of mouse pancreas at 1:500 dilution in buffer.