

Anti-FFAR3 antibody (270-320) {FITC} (STJ501016)

STJ501016

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-FFAR3 (270-320) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.58-0.59 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2865
Gene Symbol	FFAR3
Uniprot ID	FFAR3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 270-320 on human FFAR3 protein.
Immunogen Region	270-320
Specificity	This antibody will detect FFAR3. This antibody will not detect FFAR1, FFAR2, or FFAR4.
Immunogen Sequence	