

Anti-FFAR4 antibody (1-80) {FITC} (STJ501020)

STJ501020

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-FFAR4 (1-80) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.58-0.59 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	338557
Gene Symbol	FFAR4
Uniprot ID	FFAR4_HUMAN
Immunogen	N-epitope synthetic peptide taken within amino acid region 1-80 on human FFAR4 protein.
Immunogen Region	1-80
Specificity	
Immunogen Sequence	