

Anti-FGR antibody (500-529) (STJ501046)

STJ501046

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-FGR (500-529) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.61 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	2268
Gene Symbol	FGR
Uniprot ID	FGR_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 500-529 on human Fgr protein.
Immunogen Region	500-529
Specificity	
Immunogen Sequence	



Western blot with Anti-FGR Antibody (STJ501046) of T47D cell at 1:500 dilution in buffer.