

Anti-C11orf1 antibody (1-50) (STJ501058)

GENERAL INFORMATION

Product Type Primary antibodies

Short Description Rabbit polyclonal antibody anti-C11orf1 (1-50) is suitable for use in ELISA research applications.

Applications ELISA Host/Source Rabbit Reactivity Human

PRODUCT PROPERTIES

Clonality Polyclonal

Clone ID

Concentration 0.65-0.67 μg/μl Conjugation Unconjugated Purification Affinity Purified Dilution Range WB:

DB: 1:15, 000-1:20, 000 ELISA: 1:15, 000-1:20, 000

Formulation

Isotype IgG

Storage Instruction Store at-20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID 64776 Gene Symbol C11orf1 Uniprot ID CK001_HUMAN

Immunogen Synthetic peptide taken within amino acid region 1-50 on human C11orf1 protein.

Immunogen Region 1-50 Specificity Immunogen Sequence