

Anti-C11orf1 antibody (1-50) {Biotin} (STJ501061)

STJ501061

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-C11orf1 (1-50) is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.65-0.67 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: DB: 1:15, 000-1:20, 000 ELISA: 1:15, 000-1:20, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	64776
Gene Symbol	C11orf1
Uniprot ID	CK001_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-50 on human C11orf1 protein.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	