

## Anti-C11orf1 antibody (100-150) {Biotin} (STJ501063)

STJ501063

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-C11orf1 (100-150) is suitable for use in ELISA research applications.
<b>Applications</b>	ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.65-0.67 µg/µl
<b>Conjugation</b>	Biotin
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: DB: 1:15, 000-1:20, 000 ELISA: 1:15, 000-1:20, 000
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	64776
<b>Gene Symbol</b>	C11orf1
<b>Uniprot ID</b>	CK001_HUMAN
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 100-150 on human C11orf1 protein.
<b>Immunogen Region</b>	100-150
<b>Specificity</b>	
<b>Immunogen Sequence</b>	