

Anti-FTH1 antibody (Full Length) (STJ501085)

STJ501085

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-FTH1 (Full Length) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200 IHC: 1:50-1:150 ICC: 1:50-1:150
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	2495
Gene Symbol	FTH1
Uniprot ID	FRIH_HUMAN
Immunogen	Full length Ferritin Protein.
Immunogen Region	Full Length
Specificity	
Immunogen Sequence	



Western blot with Anti-FTH1 Antibody (STJ501085) of purified human Ferritin at 1:250 dilution in buffer.