

Anti-FUBP1 antibody (600-644) {Biotin} (STJ501095)

STJ501095

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-FUBP1 (600-644) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.53-0.55 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	8880
Gene Symbol	FUBP1
Uniprot ID	FUBP1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 600-644 on human FUBP1 protein.
Immunogen Region	600-644
Specificity	
Immunogen Sequence	