

Anti-F13A1 antibody (1-50) {FITC} (STJ501104)

STJ501104

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-F13A1 (1-50) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Cat, Human, Monkey, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2162
Gene Symbol	F13A1
Uniprot ID	F13A_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-50 on human factor XIII A subunit protein.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	