

Anti-F13A1 antibody (1-50) {Biotin} (STJ501107)

STJ501107

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-F13A1 (1-50) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunohistochemistry, Immunofluorescence, Immunoprecipitation and Western Blot research applications.
Applications	CM, ELISA, ICC, IHC, IF, IP, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150 IHC: 1:200 ICC: 1:200 IF: 1:200 CFM: 1:200
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	2162
Gene Symbol	F13A1
Uniprot ID	F13A_HUMAN
Immunogen	Multi-epitope synthetic peptides were used: amino acid range 1-50 and 700-732. A mixture of synthetic peptides used for FXIIIa-101AP and FXIIIa-112AP.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	