

Anti-GABRB2 antibody (350-400) {Biotin} (STJ501122)

STJ501122

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-GABRB2 (350-400) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Cat, Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2561
Gene Symbol	GABRB2
Uniprot ID	GBRB2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 350-400 on human gamma-aminobutyric acid receptor subunit beta-2 protein.
Immunogen Region	350-400
Specificity	
Immunogen Sequence	