

Anti-GAD2 antibody (50-110) {FITC} (STJ501126)

STJ501126

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-GAD2 (50-110) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.62 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	2572
Gene Symbol	GAD2
Uniprot ID	DCE2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 50-110 on human Glutamate decarboxylase 2 protein.
Immunogen Region	50-110
Specificity	
Immunogen Sequence	