

Anti-GADD45B antibody (100-160) (STJ501127)

STJ501127

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-GADD45B (100-160) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.6 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	17873
Gene Symbol	Gadd45b
Uniprot ID	GA45B_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 100-160 on mouse GADD45b protein.
Immunogen Region	100-160
Specificity	
Immunogen Sequence	



Western blot with Anti-GADD45B Antibody (STJ501127) with 1) GADD45B Positive control 2) Raw cells. Antibody at 1:500 dilution in buffer.