

Anti-GADD45B antibody (100-160) {Biotin} (STJ501128)

STJ501128

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-GADD45B (100-160) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.6 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	17873
Gene Symbol	Gadd45b
Uniprot ID	GA45B_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 100-160 on mouse GADD45b protein.
Immunogen Region	100-160
Specificity	
Immunogen Sequence	