

Anti-DDIT3 antibody (60-160) {Biotin} (STJ501131)

STJ501131

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-DDIT3 (60-160) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.57-0.68 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1649
Gene Symbol	DDIT3
Uniprot ID	DDIT3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 60-160 on human GADD153 protein
Immunogen Region	60-160
Specificity	
Immunogen Sequence	