

Anti-GADD34 antibody (380-400) (STJ501133)

STJ501133

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-GADD34 (380-400) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Bovine, Human, Monkey, Panda

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.58 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:800-1:1,000 ELISA: 1:800-1:1,000 IP: 1:200 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	23645
Gene Symbol	PPP1R15A
Uniprot ID	PR15A_HUMAN
Immunogen	Synthetic peptide corresponding to internal sequence of amino acids 380-400 of Human GADD34, covalently modified and coupled to a carrier protein.
Immunogen Region	380-400
Specificity	
Immunogen Sequence	



Immunohistochemistry of rat tissue with Anti-GADD34 primary antibody (STJ501133) at 1:100 dilution in Buffer. Section was treated with DAB and Haematoxylin stain. Magnification is at 40X.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081