

Anti-EIF2AK4 antibody (1600-1700) (STJ501136)

STJ501136

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-EIF2AK4 (1600-1700) is suitable for use in Confocal Microscopy, ELISA, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	CM, ELISA, IF, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.58 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:100 ICC: 1:100 IF: 1:100 CFM: 1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	440275
Gene Symbol	EIF2AK4
Uniprot ID	E2AK4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1600-1700 on human GCN2 protein.
Immunogen Region	1600-1700
Specificity	
Immunogen Sequence	



Immunohistochemistry of rat tissue with Anti-EIF2AK4 primary antibody (STJ501136) at 1:100 dilution in Buffer. Section was treated with DAB and Haematoxylin stain. Magnification is at 40X.