

## Anti-EIF2AK4 antibody (1600-1700) {Biotin} (STJ501140)

STJ501140

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Chicken polyclonal antibody anti-EIF2AK4 (1600-1700) is suitable for use in Confocal Microscopy, ELISA, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	CM, ELISA, IF, IHC, IP, WB
<b>Host/Source</b>	Chicken
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.58 µg/µl
<b>Conjugation</b>	Biotin
<b>Purification</b>	Affinity Purified
<b>Dilution</b>	WB: 1:500
<b>Range</b>	DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:100 ICC: 1:100 IF: 1:100 CFM: 1:100
<b>Formulation</b>	
<b>Isotype</b>	IgY
<b>Storage</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">440275</a>
<b>Gene Symbol</b>	<a href="#">EIF2AK4</a>
<b>Uniprot ID</b>	<a href="#">E2AK4_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 1600-1700 on human GCN2 protein.
<b>Immunogen Region</b>	1600-1700
<b>Specificity</b>	
<b>Immunogen Sequence</b>	