

Rabbit Anti-Goat Heavy Chain antibody {Biotin} (STJ501161)

STJ501161

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	Biotin-conjugated rabbit polyclonal anti-Goat Heavy Chain secondary antibody. For use in most research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Goat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Rabbit anti-Goat IgG (H)
Immunogen Region	
Specificity	
Immunogen Sequence	