

Anti-KCNJ9 antibody {Biotin} (STJ501170)

STJ501170

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-KCNJ9 is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.59 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3765
Gene Symbol	KCNJ9
Uniprot ID	KCNJ9_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid region on human GIRK3 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	