

Anti-KCNJ5 antibody (350-450) {FITC} (STJ501174)

STJ501174

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-KCNJ5 (350-450) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.53-0.57 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:250 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3762
Gene Symbol	KCNJ5
Uniprot ID	KCNJ5_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 350-450 on human GIRK4 protein.
Immunogen Region	350-450
Specificity	
Immunogen Sequence	