

Anti-SLC2A1 antibody {FITC} (STJ501198)

STJ501198

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Sheep polyclonal antibody anti-SLC2A1 is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	CM, ELISA, ICC, IF, IHC, IP, WB
Host/Source	Sheep
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:5, 000 DB: 1:12, 000 ELISA: 1:12, 000 IP: 1:200 IHC: 1:200 ICC: 1:200 IF: 1:200 CFM: 1:200
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	6513
Gene Symbol	SLC2A1
Uniprot ID	GTR1_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequences on GLUT1 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	