

Anti-SLC2A5 antibody {Biotin} (STJ501221)

STJ501221

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-SLC2A5 is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications. |
| Applications | ELISA, IP, WB |
| Host/Source | Rabbit |
| Reactivity | Chimpanzee, Human, Rat |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.7 µg/µl |
| Conjugation | Biotin |
| Purification | Affinity Purified |
| Dilution Range | WB: 1:1, 000 DB: 1:4, 000 ELISA: 1:4, 000 IP: 1:200 |
| Formulation | |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 6518 |
| Gene Symbol | SLC2A5 |
| Uniprot ID | GTR5_HUMAN |
| Immunogen | Synthetic peptide corresponding to unique amino acid sequence on human GLUT 5 protein. |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |