

## Anti-SLC2A9 antibody (50-570) {Biotin} (STJ501230)

STJ501230

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-SLC2A9 (50-570) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.5 µg/µl
<b>Conjugation</b>	Biotin
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500-1:750 DB: 1:4, 500 ELISA: 1:4, 500 IP: 1:200
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	56606
<b>Gene Symbol</b>	SLC2A9
<b>Uniprot ID</b>	GTR9_HUMAN
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 50-570 on Glut 9 protein.
<b>Immunogen Region</b>	50-570
<b>Specificity</b>	
<b>Immunogen Sequence</b>	