

Anti-GOLGA4 antibody (2025-2075) {Biotin} (STJ501242)

STJ501242

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-GOLGA4 (2025-2075) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2803
Gene Symbol	GOLGA4
Uniprot ID	GOGA4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 2025-2075 on human Golgin-245 protein.
Immunogen Region	2025-2075
Specificity	
Immunogen Sequence	