

Anti-GOLGA4 antibody (2025-2075) {FITC} (STJ501243)

STJ501243

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-GOLGA4 (2025-2075) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications. |
| Applications | ELISA, IHC, IP, WB |
| Host/Source | Rabbit |
| Reactivity | Human, Monkey |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.55-0.75 µg/µl |
| Conjugation | FITC |
| Purification | Affinity Purified |
| Dilution Range | WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150 |
| Formulation | |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 2803 |
| Gene Symbol | GOLGA4 |
| Uniprot ID | GOGA4_HUMAN |
| Immunogen | Synthetic peptide taken within amino acid region 2025-2075 on human Golgin-245 protein. |
| Immunogen Region | 2025-2075 |
| Specificity | |
| Immunogen Sequence | |