

Anti-GOLGA5 antibody (660-710) {FITC} (STJ501246)

STJ501246

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-GOLGA5 (660-710) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	9950
Gene Symbol	GOLGA5
Uniprot ID	GOGA5_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 660-710 on human Golgin subfamily A member 5 protein.
Immunogen Region	660-710
Specificity	
Immunogen Sequence	