

Anti-GPR15 antibody (320-360) {FITC} (STJ501252)

STJ501252

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-GPR15 (320-360) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2838
Gene Symbol	GPR15
Uniprot ID	GPR15_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 320-360 on human GPCR GPR15 protein.
Immunogen Region	320-360
Specificity	
Immunogen Sequence	