

Anti-GPCR GPR3 antibody (5-50) (STJ501253)

STJ501253

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-GPCR GPR3 (5-50) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.65 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:750 DB: 1:20, 000-1:50, 000 ELISA: 1:20, 000-1:50, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	266769
Gene Symbol	Gpr3
Uniprot ID	GPR3_RAT
Immunogen	Synthetic peptide taken within amino acid region 5-50 on rat GPCR3 protein.
Immunogen Region	5-50
Specificity	
Immunogen Sequence	



Western blot with Anti-GPCR GPR3 Antibody (STJ501253) with peptide on at 1:500 dilution in buffer.