

Anti-NPBWR1 antibody (1-50) {FITC} (STJ501261)

STJ501261

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-NPBWR1 (1-50) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10,000 IP: 1:150 IHC: 1:50-1:100
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2831
Gene Symbol	NPBWR1
Uniprot ID	NPBW1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-50 on human neuropeptides B/W receptor type 1 protein.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	