

## Anti-GPR75 antibody (25-75) {Biotin} (STJ501267)

STJ501267

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-GPR75 (25-75) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunohistochemistry, Immunofluorescence, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	CM, ELISA, ICC, IHC, IF, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.75-1.25 µg/µl
<b>Conjugation</b>	Biotin
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:250 IHC: 1:250 ICC: 1:250 IF: 1:250 CFM: 1:250
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">10936</a>
<b>Gene Symbol</b>	<a href="#">GPR75</a>
<b>Uniprot ID</b>	<a href="#">GPR75_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 25-75 on human G protein-coupled receptor 75 protein.
<b>Immunogen Region</b>	25-75
<b>Specificity</b>	
<b>Immunogen Sequence</b>	