

## Anti-GAPDH antibody (75-150) {Biotin} (STJ501278)

STJ501278

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-GAPDH (75-150) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, IHC, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.64-0.72 µg/µl
<b>Conjugation</b>	Biotin
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:100
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">14433</a>
<b>Gene Symbol</b>	<a href="#">Gapdh</a>
<b>Uniprot ID</b>	<a href="#">G3P_MOUSE</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 50-100 on mouse Glyceraldehyde-3-phosphate dehydrogenase protein.
<b>Immunogen Region</b>	75-150
<b>Specificity</b>	
<b>Immunogen Sequence</b>	