

Anti-GPS2 antibody (250-300) {Biotin} (STJ501284)

STJ501284

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-GPS2 (250-300) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2874
Gene Symbol	GPS2
Uniprot ID	GPS2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 250-300 on human G protein pathway suppressor 2 protein.
Immunogen Region	250-300
Specificity	
Immunogen Sequence	