

Anti-GZMB antibody (1-50) {Biotin} (STJ501291)

STJ501291

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-GZMB (1-50) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.52 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3002
Gene Symbol	GZMB
Uniprot ID	GRAB_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-50 on human endogenous Granzyme B protein.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	