

Anti-GAST antibody (50-100) {Biotin} (STJ501302)

STJ501302

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-GAST (50-100) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	CM, ELISA, ICC, IF, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150 IHC: 1:50-1:200 ICC: 1:50-1:200 IF: 1:50-1:200 CFM: 1:50-1:200
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	2520
Gene Symbol	GAST
Uniprot ID	GAST_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 50-100 on human Gastrin preprotein.
Immunogen Region	50-100
Specificity	
Immunogen Sequence	