

Anti-GTP cyclohydrolase 1 antibody (1-50) {Biotin} (STJ501306)

STJ501306

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-GTP cyclohydrolase 1 (1-50) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.68-0.98 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:750 DB: 1:20, 000 ELISA: 1:20, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	14528
Gene Symbol	Gch1
Uniprot ID	GCH1_MOUSE
Immunogen	N-epitope synthetic peptide taken within amino acid region 1-50 on GTPCHI protein.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	