

Anti-HIPK2 antibody (1150-1198) (STJ501343)

STJ501343

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-HIPK2 (1150-1198) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	28996
Gene Symbol	HIPK2
Uniprot ID	HIPK2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1150-1198 on human Homeodomain-interacting protein kinase 2 protein.
Immunogen Region	1150-1198
Specificity	
Immunogen Sequence	