

Anti-HIPK3 antibody (860-910) (STJ501346)

STJ501346

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-HIPK3 (860-910) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10114
Gene Symbol	HIPK3
Uniprot ID	HIPK3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 860-910 on human Homeodomain-interacting protein kinase 3 protein.
Immunogen Region	860-910
Specificity	
Immunogen Sequence	