

Anti-HIPK3 antibody (860-910) {Biotin} (STJ501347)

STJ501347

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-HIPK3 (860-910) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications. |
| Applications | ELISA, IP, WB |
| Host/Source | Rabbit |
| Reactivity | Human, Mouse, Rat |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.5-1.5 µg/µl |
| Conjugation | Biotin |
| Purification | Affinity Purified |
| Dilution Range | WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200 |
| Formulation | |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 10114 |
| Gene Symbol | HIPK3 |
| Uniprot ID | HIPK3_HUMAN |
| Immunogen | Synthetic peptide taken within amino acid region 860-910 on human Homeodomain-interacting protein kinase 3 protein. |
| Immunogen Region | 860-910 |
| Specificity | |
| Immunogen Sequence | |