

Anti-MT-RNR2 antibody (5-24) {Biotin} (STJ501372)

STJ501372

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-MT-RNR2 (5-24) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	NA
Gene Symbol	MT-RNR2
Uniprot ID	HUNIN_HUMAN
Immunogen	Synthetic peptide corresponding to amino acid region 5-24 on Humanin protein.
Immunogen Region	5-24
Specificity	Will detect all isoforms except Humanin 9.
Immunogen Sequence	