

Anti-B19 antibody (650-714) (STJ501374)

STJ501374

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-B19 (650-714) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200 IHC: 1:50-1:100
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	26559
Gene Symbol	Hunk
Uniprot ID	HUNK_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 650-714 on mouse Hormonally Up-regulated Neu tumor-associated Kinase.
Immunogen Region	650-714
Specificity	
Immunogen Sequence	