

Anti-IDO1 antibody (100-150) {Biotin} (STJ501390)

STJ501390

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-IDO1 (100-150) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3620
Gene Symbol	IDO1
Uniprot ID	I23O1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 100-150 on human Indoleamine 2, 3-dioxygenase 1 protein.
Immunogen Region	100-150
Specificity	
Immunogen Sequence	