

## Anti-IFI16 antibody (140-190) (STJ501395)

STJ501395

### GENERAL INFORMATION


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-IFI16 (140-190) is suitable for use in ELISA and Western Blot research applications.
<b>Applications</b>	ELISA, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.58-0.62 µg/µl
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	3428
<b>Gene Symbol</b>	IFI16
<b>Uniprot ID</b>	IF16_HUMAN
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 140-190 on human IFI16 protein
<b>Immunogen Region</b>	140-190
<b>Specificity</b>	
<b>Immunogen Sequence</b>	



Western blot with Anti-IFI16 Antibody (STJ501395) of THP cell at 1:500 dilution in buffer.