

Anti-IFI16 antibody (140-190) {Biotin} (STJ501396)

STJ501396

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-IFI16 (140-190) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.58-0.62 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3428
Gene Symbol	IFI16
Uniprot ID	IF16_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 140-190 on human IFI16 protein
Immunogen Region	140-190
Specificity	
Immunogen Sequence	