

Anti-IFITM3 antibody (1-50) {Biotin} (STJ501405)

STJ501405

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-IFITM3 (1-50) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:150 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10410
Gene Symbol	IFITM3
Uniprot ID	IFM3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-50 on human Interferon-induced transmembrane protein 3.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	